



Elemental Analysis: NCS characterization of fertilizers

Authors

Dr. Liliana Krotz and
Dr. Guido Giazzi
Thermo Fisher Scientific,
Milan, Italy

Keywords

Automated Analysis, Combustion,
Fertilizers, NCS, Unattended
Analysis

Goal

This application notes reports data of different solid and liquid fertilizer samples to show the accuracy of the FlashSmart Elemental Analyzer for samples with large range of element concentrations.

Introduction

Fertilizers are chemically characterized for research and quality control purposes in agronomy.

Nitrogen and carbon determination in soils enables the evaluation of organic matter and the calculation of the needed amount of fertilizer since characterizing nitrogen and carbon provides information on nutritional elements of the organic matter. When a lack of sulfur is determined, growth conditions of vegetables is negatively affected, and consequently the quality of the protein that they produce. For this reason the demand for the determination of nitrogen, carbon and sulfur in fertilizers and in agricultural water resulting from the use of fertilizer has increased. It is therefore very important to have an accurate and automated analytical technique, to perform fast analysis with excellent reproducibility.

The Thermo Scientific™ FlashSmart™ Elemental Analyzer (Figure 1) based on the dynamic combustion method (modified Dumas method), provides automated nitrogen, carbon and sulfur determination without the use of hazardous chemicals. The FlashSmart™ EA enables to run determinations at high and low levels of concentrations for solid and liquid samples in one single system.

This note presents data on NCS determination of fertilizers to show the performance of the FlashSmart Elemental Analyzer.

Methods

The FlashSmart EA operates according to the dynamic flash combustion of the sample. Samples are weighed in tin containers and introduced into the combustion reactor via the Thermo Scientific™ MAS Plus Autosampler or they can be directly injected by a syringe via the Thermo Scientific™ AS 1310 Liquid Autosampler, in both cases with oxygen.

Solid samples were homogenized by a ball mill while liquid samples were analyzed without pre-treatment.

For nitrogen only determination, after combustion, the resultant gases are carried by a helium flow to a second reactor filled with copper, then swept through CO₂ and H₂O traps, a GC column and finally detected by a Thermal Conductivity Detector (TCD). Total run time less than five minutes (Figure 2).

For simultaneous NCS or for sulfur only determination, after combustion, the gases are carried by a helium flow to a layer filled with copper, then swept through a water trap, a GC column that provides the separation of the combustion gases, and finally, detected by a Thermal Conductivity Detector. Total run time is 5 – 10 min. (Figure 3).

A complete report is automatically generated by the Thermo Scientific™ EagerSmart™ Data Handling Software and displayed at the end of the analysis.

Results

This paper presents data of different solid and liquid fertilizer samples analyzed several times. Liquid samples were homogenized by a ball mill while liquid samples were analyzed without pre-treatment.

Table 1 shows the nitrogen determination of solid fertilizer samples. Instrument calibration was performed with approximately 100 mg of urea (46.65 N%). Urea was analyzed also as an unknown sample to check the accuracy and precision of the system. The weight used for fertilizers was 4 – 5 mg for ammonium nitrate, 100 - 120 mg for potassium nitrate and 250 - 270 mg for the mix.



Figure 1. Thermo Scientific FlashSmart Elemental Analyzer.




Figure 2. Nitrogen determination.




Figure 3. NCS or sulfur determination.

Table 1. Nitrogen determination in solid fertilizers.

Sample	N%	RSD%
Urea	46.59 - 46.56 - 46.56	0.268
Thermo Scientific Standard 46.67 N%		
Potassium nitrate KNO ₃	13.80 - 13.82 - 13.94	0.561
Ammonium nitrate NH ₄ NO ₃	34.93 - 34.93 - 34.93	0.009
Urea and ammonium sulphate (NH ₄) ₂ SO ₄ mix	34.56 - 34.71 - 34.78	0.324
Fertilizer mix A	17.55 - 17.51 - 17.55	0.147
Fertilizer mix B	15.44 - 15.42 - 15.57	0.524
Fertilizer mix C	39.055 - 39.131 - 38.895	0.309
Fertilizer mix D	38.765 - 38.821 - 38.729	0.121

Table 2 shows the nitrogen determination of liquid fertilizer samples. Samples were analyzed by liquid

injection at 100 µl. Instrument calibration was performed with 100 ul of urea water solutions (0.3 and 7 N%).

Table 2. Nitrogen determination in liquid fertilizers.

Sample	N%	RSD%	Sample	N%	RSD%
1	0.072	1.519	6	2.789	1.075
	0.073			2.849	
	0.072			2.830	
2	0.220	0.607	7	7.133	0.324
	0.218			7.103	
	0.220			7.148	
3	0.310	1.103	8	10.197	0.193
	0.316			10.160	
	0.316			10.167	

Table 3 shows the nitrogen and carbon determination of solid fertilizers. To calibrate the system, 4 - 5 mg of

aspartic acid (10.52 N%, 36.09 C%) was used and samples were weighed at 5 - 6 mg.

Table 3. Nitrogen and carbon determination in solid fertilizers.

Sample Name	Nitrogen		Carbon	
	%	RSD%	%	RSD%
1	3.876	0.274	29.837	1.250
	3.869		30.476	
	3.891		29.813	
2	17.873	0.215	0.231	0.910
	17.807		0.227	
	17.874		0.228	
3	33.466	0.125	0.038	1.860
	33.505		0.037	
	33.550		0.039	

Table 4 shows the nitrogen and carbon determination of liquid fertilizers. To calibrate the system 4 - 5 mg of aspartic acid (10.52 N%, 36.09 C%) was used and

samples were weighed at 10 - 15 mg adsorbed on approximately 10 mg of Chromosorb® WAW 30-60 mesh.

Table 4. Nitrogen and carbon determination in liquid fertilizers.

Sample	N%	RSD%	C%	RSD%
1	0.7881	0.2412	0.4552	1.8852
	0.7752		0.4675	
2	1.6335	0.0260	0.0133	6.1005
	1.6329		0.0122	
3	4.0065	0.9506	0.0416	2.0696
	3.9530		0.0404	
4	4.0576	0.8256	0.3475	1.5497
	4.0105		0.3552	
5	0.6753	1.3483	0.4791	0.1180
	0.6883		0.4799	

Table 5 shows the NCS data of solid fertilizer samples. The calibration was performed with 2 - 3 mg of BBOT*, Nicotinamide, Cystine and Urea, and the sample weight used was 3 - 4 mg.

* BBOT (6.51 N%, 72.53 C%, 7.44 S%), Nicotinamide (22.94 N%, 59.01 C%), Cystine (11.66 N%, 29.99 C%, 26.69 S%) and Urea (46.65 N%, 20 C%)

Table 5. Simultaneous NCS determination in solid fertilizers.

Sample	N%	RSD%	C%	RSD%	S%	RSD%
1	10.538	0.230	4.583	1.821	0.218	1.862
	10.582		4.523		0.222	
	10.578		4.688		0.214	
2	24.613	0.247	0.221	0.686	1.232	0.126
	24.732		0.223		1.234	
	24.698		0.220		1.231	
3	25.500	0.094	0.046	3.901	0.085	0.713
	25.538		0.042		0.084	
	25.546		0.044		0.085	
4	20.632	0.305	0.943	0.411	23.939	0.115
	20.755		0.936		23.986	
	20.669		0.939		23.987	

Table 6 shows the sulfur data of solid samples. The calibration was performed with 2 - 3 mg of BBOT* (7.44 S%) and the sample weight used was 2 - 4 mg.

* BBOT: 2,5-Bis (5-tert-butyl-benzoxazol-2-yl) thiophene

Table 6. Sulfur determination in solid fertilizers.

Fertilizer 1		Fertilizer 2		Fertilizer 3	
S%	RSD%	S%	RSD%	S%	RSD%
2.740	0.735	3.616	0.857	4.873	0.638
2.698		3.614		4.849	
2.722		3.664		4.860	
2.749		3.623		4.810	
2.738		3.681		4.803	

Conclusions

All data were obtained with a good reproducibility and no matrix effect was observed when changing from solid to liquid samples.

The FlashSmart EA enables to perform simultaneous NCS determination in a single run, and consequently sulfur, with no hardware modification required.

Nitrogen only determination can be performed by increasing the sample weight and changing the configuration.

The FlashSmart EA allows to obtain automated analysis with lower cost per sample, while obtaining excellent accuracy and reproducible data.

Find out more at thermofisher.com/OEA

ThermoFisher
SCIENTIFIC